

AMP HANDRIM

Mark For: _____
 Acct #: _____
 Bill To Name: _____
 Billing Address: _____
 City, State: _____
 Confirmation Email: _____
 Contact : _____

QUOTE ☐ ORDER ☐ #: _____
 Date: _____
 PO#: _____
 Ship to: _____
 Address: _____
 City, State: _____
 Zip: _____



END USER WHEELCHAIR SPECIFICATIONS

Wheelchair Type: Folding ☐ Rigid ☐ Brand: _____ Model: _____
 1. Front Seat Height in inches  _____
 2. Rear Seat Height in inches  _____
 3. Inside Tire Distance in inches  _____

¹ Wheelchair front seat height must be a minimum 1/2" greater than rear; to maintain same incline when mounted to AMP add an additional 1/2" in front height of wheelchair as docking increases rear seat height 1/2". ² To avoid wheelchair interference with the AMP Module the rear lower frame must have 8.5" ground clearance. ³ At bottom of chair inside to inside distance between tires.

AMP HANDRIM & OPTIONS

- ☐ AMP E0986 AMP Handrim Power Module w/ Battery & Charger
 (Must choose one from the size matrix to the right:)
☐ 380AHRS AMP Small Base 14" - 16" Wide
☐ 380AHM AMP Medium Base 16" - 18" Wide
☐ 380AHL AMP Large Base 18" - 20" Wide
☐ 180WEXT Width Extension Kit (Adds 2" to Wheel Width Fit)

Joystick Input

- ☐ 380JIO Joystick Input
 Joystick Mount Side (Must choose one if ordering Joystick option)
☐ 380JIR Joystick Input Right
☐ 380JIL Joystick Input Left
 Swing Away/Attendant Control (Must choose at least one if ordering Joystick option)
☐ 380SAD Swing Away Mount Dual Post (Desk Pad Style)
☐ 380SAT Swing Away Mount Tubular

Joystick Handles

- ☒ 380SH Standard Handle (When ordering Joystick input)
☐ 380BPU Bodypoint® 3" U-Shaped Handle
☐ 380BPGM Bodypoint® Grooved Mushroom Handle

Battery/ Charger

- ☒ 380LPSTC LiFePO4 Battery/ Charger (6.5 Hours)
☐ 380USC Upgrade Super Charger (2.5 hours)

Battery Options

- ☐ 380BULP Back Up LiFePO4 Battery
☐ 380BUSEC Spare Battery Charger (6.5 Hours)
☐ 380SSC Spare Super Charger (2.5 Hours)

Sub Totals

AMP SIZE MATRIX

NOTE: Choose AMP size based on chair type, along with camber, which determines the AMP size recommended. Each size is adjustable in width, but limited based on camber or overall wheel base width.

CHECK ONE BELOW

FOLDING/RIGID 0/3 CAMBER		WHEEL WIDTH RANGE
CHAIR SEAT WIDTH	AMP SIZE (growth)	INSIDE TO INSIDE OF TIRES
14"	☐ SMALL (14"-16")	17 1/4" - 19 1/2"
15"	☐ SMALL (14"-16")	
16"	☐ SMALL (14"-16")	
16"	☐ MEDIUM (16"-18")	19" - 21 1/4"
17"	☐ MEDIUM (16"-18")	
18"	☐ MEDIUM (16"-18")	
18"	☐ LARGE (18"-20")	20 1/2" - 23"
19"	☐ LARGE (18"-20")	
20"	☐ LARGE (No Growth)	

FOLDING/RIGID UP TO 5 CAMBER		WHEEL WIDTH RANGE
CHAIR SEAT WIDTH	AMP SIZE (growth)	INSIDE TO INSIDE OF TIRES
14"	☐ SMALL (14"-15")	17 1/4" - 19 1/2"
15"	☐ MEDIUM (15"-17")	
16"	☐ MEDIUM (15"-17")	
17"	☐ LARGE (17"-19")	19" - 21 1/4"
18"	☐ LARGE (17"-19")	
19"	☐ LARGE (No Growth)	
20"	N/A	



For quotes and orders EMAIL to:
sales@paragonmobility.com.au

Method Mobility VA Registration: UVRSAU9M5AE7/CAGE #9JCF9

www.methodmobility.com
www.paragonmobility.com.au

Shipping
 Total w/Delivery & Shipping
 Total with Fees (Pay this amount)

